

www.phantompower.com.au

tech support ☎ 1300 667 224

Jeannie Crusher Operation



- A new bag **MUST** be used for each use of the machine.
- Insert pills into bag (recommended maximum of 6 pills per bag) with Jeannie Logo facing upwards and arrows towards machine.
- Insert bag into Jeannie in the direction of the arrows. Operation will start automatically.
- **Tip:** Pull bag through slowly to avoid jamming and damage to pill bag.
- **Tip:** View our training videos for more information at www.phantompower.com.au



Insert bag. Arrows on bag must match arrows on machine.



Take other side with opposite hand.



Pull gently through to other side. should take 3 to 5 seconds.

www.phantompower.com.au

☎ 1300 667 224



www.phantompower.com.au
tech support  **1300 667 224**



Charging

- It is recommended to leave the Jeannie Crusher on charge when not in use.
- The Jeannie Crusher should be recharged every day. This can be done over night to ensure a full charge for the next day.
- View our helpful training video on the correct charging technique at www.phantompower.com.au

Jeannie Crusher Indicator Light

The Jeannie Crusher has an indicator light to advise the status of it's operation.

- Indicator Light on constantly - Crusher is in operation
- Indicator Light single flash - Crusher needs recharging
- Indicator Light double flash - Crusher is jammed
- View our training videos on troubleshooting at www.phantompower.com.au

Maintenance

- It is recommended to keep the bag entry point free from dust and liquids by cleaning with a clean damp cloth.
- Return to Phantom Power Products for service **annually**. Call **1300 667 224**

Trouble Shooting

- If you are having problems with the operations your Jeannie Pill Crusher, please follow our 3 step troubleshooting guidelines below;
 1. View our operation and troubleshooting guides online at www.phantompower.com.au
 2. Call our technical support hotline on **1300 667 224**
 3. If you are still experiencing difficulties, return the machine to us at Phantom Power Products for repair. Please include proof of purchase, and contact details for return.

www.phantompower.com.au



1300 667 224